## Manual



Article	Content	
SL-9802-smp	1 mL, 100 rxn á 25µL	
SL-9802-5ML	4x 1.25 mL, 400 rxn á 25μL	
SL-9802-10ML	8x 1.25 mL, 800 rxn á 25μL	
SL-9802-20ML	16 x 1.25 mL, 1600 rxn á 25μL	

# primaQUANT PROBE


# 2x qPCR Probe Master Mix - no ROX

#### **DESCRIPTION**

Our **primaQUANT** *PROBE* **2x qPCR Master Mix** is an optimized ready-to-use mixture for probe-based assays such as Taqman®, Beacons and MGBs. It contains a modified HotStart DNA Polymerase, as well as dNTPS and MgCl<sub>2</sub> - combined in an optimized buffer system that provides fast kinetics and target amplification even for difficult templates.

The **primaQUANT 2x Master Mix** contains all components - you just need to add primers and template DNA/cDNA.

The Master Mix can be used not only for expression analysis but also for genotyping, copynumber analysis and all sorts of probe-based quantitative PCR.



- **primaQUANT** *PROBE* is also offered as a **blue**-colored mix for better handling.
- Some qPCR cyclers require ROX **primaQUANT PROBE** is also available with low or high concentrations of ROX.



#### STANDARD PROTOCOL



- After thawing, please **invert the Master Mix tube 6-8 times.**
- **Do not vortex** the Master Mix to prevent damage of the enzyme.

# (i) NOTE

- > Cycling conditions highly depend on the primer, probe, amplicon and input material and thus might need adjustments.
- However, standard cycling conditions can be applied for the majority of qPCR assays out-of-the box.

## 3-Step Protocol

Step	Time	Temperature	
Initial Denaturation	1-3 minutes 92 °C - 95 °C		
25- 40 cycles			
Denaturation	5 - 10 seconds 92 °C - 95 °C		
Annealing	5 seconds 60 °C - <b>depending on primer</b>		
Extension	<b>20 - 30 seconds</b> 72 °C		

## 2-Step Protocol

Step	Step Time		
Initial Denaturation	1-3 minutes 92 °C - 95 °C		
25- 40 cycles			
Denaturation	5 seconds 92 °C - 95 °C		
Annealing/Extension Combined	20 - 40 seconds	60 °C - depending on primer	



#### **ULTRA-FAST PROTOCOL**

- (i) BEFORE YOU START
  - After thawing, please **invert the Master Mix tube 6-8 times.**
  - **Do not vortex** the Master Mix to prevent damage of the enzyme.
- (i) NOTE
  - Ultra-fast Cycling conditions highly depend on the ramping rate of your qPCR cycler, primer, probe, amplicon and input material and thus might need adjustments.
  - Ultra-fast cycling conditions can be applied for the majority of qPCR assays out-of-the box, provided that your primer/probe sets do not show unspecific binding.

# 3-Step Protocol

Step	Time	Temperature	
Initial Denaturation	60 seconds	92 °C - 95 °C	
25- 40 cycles			
Denaturation	1-5 seconds	92 °C - 95 °C	
Annealing	1-5 seconds	60 °C - depending on primer	
Extension	5 seconds 72 °C		

## 2-Step Protocol

Step	Time	Temperature	
Initial Denaturation	60 seconds	92 ℃ - 95 ℃	
25- 40 cycles			
Denaturation	1 second	92 ℃ - 95 ℃	
Annealing/Extension Combined	5-10 seconds 60 °C - depending on p		



### RECOMMENDED REACTION MIXTURE PER WELL

Components	25 μL Reaction	<b>10</b> μL Reaction	Final Concentration
2x primaQUANT Master Mix	12.5 µL	5 μL	1x
Forward Primer	variable (e.g. 2 µL)	variable (e.g. 1 µL)	100 - 400 nM
Reverse Primer	variable (e.g. 2 µL)	variable (e.g. 1 µL)	100 - 400 nM
Probe	variable (e.g. 2 µL)	variable (e.g. 1 µL)	200 - 600 nM
Template DNA	variable	variable	0.1 - 10 ng/reaction
Sterile Water	adjust to 25 μL	adjust to 10 μL	





For maximum efficiency and specificity annealing temperatures as well as extension time, primer/probe concentration and template concentration need to be optimized.

### **CALCULATOR TOOL**



Please feel free to download our Excel sheet calculator to calculate the necessary volumes:

http://calculator.steinbrenner-laborsysteme.de.



https://www.steinbrenner-laborsysteme.de



## HOW TO VALIDATE A QPCR SETUP



You can find additional information on how to validate and set up a qPCR in our qPCR Knowledge Center.

# **Required Controls**

## > DNA Dilution Series

A DNA dilution series is used to validate the dynamic range, find the optimal DNA input amount and estimate the overall PCR efficiency.

### > Primer Dilution Series

High primer amounts can result in unspecific primer binding that limit the fidelity of your qPCR.


## No Template Control (NTC)

A control in which all components except the template are added - this control is used as a negative control and should not show amplification.

## No Reverse Transcription (NoRT)

For this control, reverse transcription is performed without the reverse transcriptase. It helps to identify unwanted amplification from left-over gDNA or unwanted templates.

## **Recommended Validation Layout**





#### **APPLICATIONS**

Probe-based quantitative PCR

TaqMan® Probes
Any Dual-Labeled Hydrolysis Probe
Molecular Beacons
Scorpion Probes

- > DNA Genotyping
- DNA SNP Analysis
- RNA and miRNA Expression
- Multiplexing (up to 4 colors)
- Transcript Variant Analysis

## **QUALITY CONTROL**

Our **primaQUANT** *PROBE* **2x qPCR Master Mix** undergoes stringent quality controls. Each lot is tested in a probe-based qPCR with cDNA and lambda DNA input.

Enzyme purity and homogenity of > 98 % is validated using a Bioanalyzer SDS protein electrophoresis.

All Master Mixes are free of detectable endonuclease- & exonuclease activity:

- Incubation of 1 µg of plasmid DNA with 5 U for 4h at 37 °C and 72 °C
- Incubation of 1 μg of a DNA size standard with 5 U for 4h at 37 °C and 72 °C

#### **FURTHER INFORMATION**

For more information, please visit our website

https://www.steinbrenner-laborsysteme.de





Steinbrenner Laborsysteme GmbH In der Au 17 · 69257 Wiesenbach · Germany Fon +49 (0)6223 8612-47 Fax +49 (0)6223 8612-48 mail@steinbrenner-laborsysteme.de www.steinbrenner-laborsysteme.de