



## Taq DNA Polymerase with 10x PCR Buffer

Article	Content
SL-9601/02-smp	200 U, 200 rxn × 20 µl
SL-9601/02-500	500 U, 500 rxn × 20 µl
SL-9601/02-2500	5 × 500 U, 2500 rxn × 20 µl
SL-9601/02-10000	20 × 500 U, 10000 rxn × 20 µl



**Long-Term Storage**  
at -20°C in the dark

**Short-Term Storage**  
at 4°C in the dark

### DESCRIPTION

Our **primaAMP Taq DNA polymerase** is a recombinant, thermostable Taq from *Thermus aquaticus*. It possesses a 5' to 3' polymerase activity as well as a 5'-flap endonuclease activity.

The **primaAMP Taq DNA polymerase** can be used to amplify DNA fragments up to a length of 5 kb. Moreover, it generates A (adenine) overhangs at the 3' end, which can be used for TA-cloning.

The optimized 10x PCR buffer contains potassium chloride as well as ammonium sulfate and allows the amplification of difficult templates (e.g. GC-rich).

The polymerase is available with **two different 10x PCR buffers**. The **red 10x PCR Buffer includes a DNA loading dye** and can be used to directly load the PCR sample onto agarose gels after cycling.



### DID YOU KNOW?

- **primaAMP** is also available as a ready-to-use 2x PCR Master Mix (SL-9611/SL-9612).
- For HotStart polymerase, please order SL-9701/SL-9702.



## Recommended Reaction Mixture



### BEFORE YOU START

- After thawing, please **invert the component tubes 6-8 times**.
- **Do not vortex** the reaction mixture to avoid damaging the enzyme.

Component	20 $\mu$ l Reaction	10 $\mu$ l Reaction	Final Concentration
10x PCR Buffer (red or white)	2 $\mu$ l	1 $\mu$ l	1x
primaAMP Polymerase	0.2 $\mu$ l	0.1 $\mu$ l	1.25 U
Forward Primer	variable (e.g. 2 $\mu$ l)	variable (e.g. 1 $\mu$ l)	100 - 400 nM
Reverse Primer	variable (e.g. 2 $\mu$ l)	variable (e.g. 1 $\mu$ l)	100 - 400 nM
dNTPs (dATP, dCTP, dGTP, dTTP)	variable	variable	200 $\mu$ M each
Template DNA	variable	variable	0.01 - 10 ng per reaction
Sterile Water	adjust to 20 $\mu$ l	adjust to 10 $\mu$ l	-



## Suggested Cycling Conditions

Step	Time	Temperature
Initial Denaturation	3 minutes	92°C - 95°C
<b>25 - 35 cycles</b>		
Denaturation	5 - 10 seconds	92°C - 95°C
Annealing	5 - 10 seconds	55°C - 68°C <b>depends on primer</b>
Extension	5 - 30 seconds per 1 kb amplicon length	72°C



### NOTE

- The optimal annealing temperature is usually 2°C - 5°C below the primer melting temperature.
- Recommended elongation time is 5 - 30 seconds per 1 kb of amplicon length. For more complicated templates, we suggest 45 seconds for elongation.
- For maximum yield and specificity, we recommend to optimize annealing temperatures, annealing time, extension time, and the number of cycles for each template and primer pair.



## Contents

	Components	Description
SL-9601	SL-9601	Taq DNA Polymerase
	SL-9001	10x PCR Buffer, colourless
SL-9602	SL-9601	Taq DNA Polymerase
	SL-9002	10x PCR Buffer, red; ready-to-load

## Further Information

For more information, please visit our website: [www.steinbrenner.de](http://www.steinbrenner.de)



**Steinbrenner**  
Laborsysteme GmbH

In der Au 17 | 69257 Wiesench

+49 (0) 6223 / 96 73 00

mail@steinbrenner.de



## Further Products

Products that may also be interesting to you



- Direct-PCR without DNA extraction
- From sample to PCR in 15 minutes
- For cell culture, tissue, plants, mouse tails/ear, meat



- High-fidelity / proofreading PCR
- For NGS and cloning

