

Article	Reactions
SL-9130	1 ml, 100 rxn x 20µl
SL-9130-5ML	5x 1 ml, 500 rxn x 20µl
SL-9130-10ML	10x 1 ml, 1000 rxn x 20µl
SL-9130-20ML	20x 1 ml, 2000 rxn x 20µl
SL-9130-50ML	50x 1 ml, 5000 rxn x 20µl

Storage Conditions

Long-Term Storage
Master Mix at -20 °C

Short-Term Storage
Master Mix at 4 °C



Ultra-Fast Direct-PCR

DESCRIPTION

The **primaDIRECT** 2x PCR Master Mix enables direct PCR for fast and reliable genotyping and other standard PCR applications.

In combination with our **primaDIRECT** High Performance Lysis Buffer (SL-9102), **samples of various origin** (tissue, cells, mouse tails, plant pieces) are lysed within 15 minutes using our proprietary buffer **without the need of column- or magnetic-based extraction procedures**.

This 2x PCR Master Mix does not contain loading dye.

FURTHER INFORMATION

For more information, please visit our website

<https://www.steinbrenner-laborsysteme.de>



Steinbrenner
Laborsysteme GmbH

In der Au 17
69257 Wiesenbach

+49 (0) 6223 / 96 73 00
mail@steinbrenner.de
www.steinbrenner.de

SAMPLE LYSIS (NOT INLCUDED) & PCR WORKFLOW

BEFORE YOU START

- > The High Performance Lysis Buffer is not included and can be obtained seperately (SL-9102-20ML).
- > **If you use your own lysis buffer system, please continue with step 2.**

STEP 1: SAMPLE LYSIS

1. Fill 50 - 100 mg of your sample in to a 1.5 mL reaction tube
 1. If using cells, fill 10^3 - 10^6 cells in to a tube, spin it down at max speed and remove the supernatant medium
2. Add 50-200 μ l RTL buffer
3. Incubate at 56°C for 15-30 minutes
4. Incubate at 95 °C for 3-20 minutes to stop the enzymatic reactions
5. **Optional:** Spin down your sample and transfer the supernatant to a fresh 1.5 mL tube - recommended for long-term storage
6. Place your sample on ice while preparing the PCR reaction
7. Add 2 μ l of your lysed sample to the PCR reaction mix.

LYSIS OF SAMPLE MATERIAL

- > **The lysis will not completely degrade the sample - which is absolutely fine and does not have any influence on the PCR.**
- > You can lyse various types of sample material such as cells, animal tissue, plant leaves, mouse ears or mouse tails.
- > In general, use as little sample material for lysis as possible to reduce the amount of PCR inhibitors present in the reaction mixture.

STEP 2: PIPETTE PCR REACTION MIXTURE

Components	20 µl Reaction	Final Concentration
2x primaDIRECT Master Mix	10 µl	1x
Forward Primer	variable (e.g. 2 µl)	200 - 400 nM
Reverse Primer	variable (e.g. 2 µl)	200 - 400 nM
Lysed Sample	1-2 µl	-
Sterile Water	adjust to 20 µl	

STEP 3: RUN PCR

Step	Time	Temperature
Initial Denaturation	2 minutes	95 °C
35 - 40 cycles		
Denaturation	5 seconds	95 °C
Annealing	10 seconds	55 °C - 68 °C
Extension	variable, depends on fragment length 5 - 20 seconds per 1 kb fragment length	72 °C



Please see pages 5 & 6 for faster PCR protocols

STEP 4: VISUALIZE ON GEL

Transfer 20 µl to the agarose gel - **loading dye is required.**




NOTE



This **primaDIRECT** 2x Master Mix does not contain a loading dye. Additional loading dye is required for gel visualization.

GRAPHICAL WORKFLOW

LYSE + INACTIVATE
PCR



STANDARD GENOTYPING PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Cycling conditions highly depend on the primer, probe, amplicon and input material and thus might need adjustments.
- > However, standard cycling conditions can be applied for the majority of qPCR assays out-of-the box.

3-Step Protocol

Step	Time	Temperature
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 - 10 seconds	92 °C - 95 °C
Annealing	5 seconds	60 °C - depending on primer
Extension	20 - 30 seconds	72 °C

2-Step Protocol

Step	Time	Temperature
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 seconds	92 °C - 95 °C
Annealing/Extension Combined	20 - 40 seconds	60 °C - depending on primer

ULTRA-FAST GENOTYPING PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Ultra-fast Cycling conditions **highly depend on the ramping rate of your PCR cycler, primer, amplicon and input material** and thus might need adjustments.
- > Ultra-fast cycling conditions can be applied for the majority of PCR assays out-of-the box, provided that your primer/probe sets do not show unspecific binding.

3-Step Protocol

Step	Time	Temperature
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1-5 seconds	92 °C - 95 °C
Annealing	1-5 seconds	60 °C - depending on primer
Extension	5 seconds	72 °C

2-Step Protocol

Step	Time	Temperature
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1 second	92 °C - 95 °C
Annealing/Extension Combined	5-10 seconds	60 °C - depending on primer