



primaQUANT 1-Step RT-qPCR Master Mix for Probes

| Article | Content |
|--------------|-------------------------|
| SL-9560_smp | 0.5 ml, 50 rxn × 20 µL |
| SL-9560-1ML | 1 ml, 100 rxn × 20 µL |
| SL-9560-5ML | 5 ml, 500 rxn × 20 µL |
| SL-9560-10ML | 10 ml, 1000 rxn × 20 µL |
| SL-9560-20ML | 20 ml, 2000 rxn × 20 µL |



Long-Term Storage
at -20°C in the dark

Short-Term Storage
at 4°C in the dark

DESCRIPTION

Our **primaQUANT 1STEP PROBE RT-qPCR Master Mix** for probes is an optimized ready-to-use Master Mix for probe-based assays such as Taqman®, Beacons, MGB or Mediator Probes and can be directly used with both DNA and RNA as a starting material.

It contains a modified and proprietary HotStart DNA Polymerase and Reverse Transcriptase, as well as dNTPS, MgCl₂ and other components in optimized concentrations. It provides fast kinetics, a RNA sensitivity of < 10 fg, target amplification even for difficult templates and multiplexing of more than 6 targets.

The **primaQUANT 1STEP PROBE RT-qPCR 2x qPCR Master Mix** contains all components - you just need to add primers and template DNA/cDNA or RNA.



DID YOU KNOW?

For qPCR cyclers requiring ROX as a reference dye, **primaQUANT 1STEP PROBE Master Mix** is also available with ROX (SL-9560R).

Recommended Reaction Mixture per Well



BEFORE YOU START

- After thawing, please **invert the Master Mix tube 6-8 times**.
- **Do not vortex** the Master Mix to prevent damage of the enzyme.

| Component | Stock Concentration | 20 μ l Reaction | 10 μ l Reaction | Final Concentration |
|--|---------------------|----------------------|----------------------|-----------------------------------|
| 2x primaQUANT 1STEP PROBE RT-qPCR Master Mix | 2x | 10 μ l | 5 μ l | 1x |
| Forward Primer | 4 μ M | 1 μ l | 0.5 μ l | 200 nM (100 - 400 nM recommended) |
| Reverse Primer | 4 μ M | 1 μ l | 0.5 μ l | 200 nM (100 - 400 nM recommended) |
| Probe | 8 μ M | 1 μ l | 0.5 μ l | 400 nM (200 - 600 nM recommended) |
| Template RNA | - | variable | variable | 0.1 - 100 ng / Reaction |
| Steriles Wasser | - | adjust to 20 μ l | adjust to 10 μ l | - |



NOTE

For maximum efficiency and specificity, adjustments of annealing temperature as well as extension time, primer/probe concentration and template concentration may be needed.

CALCULATOR TOOL

Please feel free to download our Excel sheet calculator to calculate the necessary volumes:
calculator.steinbrenner-laborsysteme.de



qPCR
 KnowledgeCenter

Standard Protocol



NOTE

- For the majority of RT-qPCR assays, standard cycling conditions can be applied.
- However, cycling conditions strongly depend on the primer, probe, amplicon and input material and thus some of these factors might need adjustments.

ONE-STEP RT-QPCR WITH ADDITIONAL ANNEALING

| | | Time | Temperature |
|-----------------------|--------------|----------------|-----------------------------|
| Reverse Transcription | | 10 minutes | 50°C |
| Initial Denaturation | | 1 - 3 minutes | 92°C - 95°C |
| 25 - 40 cycles | Denaturation | 5 seconds | 92°C - 95°C |
| | Annealing | 5 seconds | 60°C depending on primer |
| | Extension | 5 - 10 seconds | 72°C |

ONE-STEP RT-QPCR WITHOUT ADDITIONAL ANNEALING

| | | Time | Temperature |
|-----------------------|-----------------------------------|----------------|-----------------------------|
| Reverse Transcription | | 10 minutes | 50°C |
| Initial Denaturation | | 1 - 3 minutes | 92°C - 95°C |
| 25 - 40 cycles | Denaturation | 5 seconds | 92°C - 95°C |
| | Annealing / Extension combined | 5 - 20 seconds | 60°C depending on primer |

Ultra-fast Protocol



NOTE

- Ultra-fast cycling conditions highly depend on the ramping rate of your qPCR cycler, primer, probe, amplicon and input material and thus might need adjustments.
- Ultra-fast cycling conditions can be applied for the majority of qPCR assays out-of-the box, provided that your primer/probe sets do not show unspecific binding.

ONE-STEP RT-QPCR PROTOCOL WITH ADDITIONAL ANNEALING

| | | Time | Temperature |
|-----------------------|--------------|----------------|-----------------------------|
| Reverse Transcription | | 5 - 10 minutes | 50°C |
| Initial Denaturation | | 1 minute | 92°C - 95°C |
| 25 - 40 cycles | Denaturation | 1 second | 92°C - 95°C |
| | Annealing | 1 - 5 seconds | 60°C depending on primer |
| | Extension | 1 second | 72°C |

ONE-STEP RT-QPCR PROTOCOL WITHOUT ADDITIONAL ANNEALING

| | | Time | Temperature |
|-----------------------|-----------------------------------|----------------|-----------------------------|
| Reverse Transcription | | 5 - 10 minutes | 50°C |
| Initial Denaturation | | 1 minute | 92°C - 95°C |
| 25 - 40 cycles | Denaturation | 1 second | 92°C - 95°C |
| | Annealing / Extension combined | 1 - 5 seconds | 60°C depending on primer |

Applications

Probe-based quantitative PCR with RNA Input

Probe-based quantitative PCR with cDNA/DNA Input

Multiplex RT-qPCR Assays for up to 8 colors

QUALITY CONTROL

Our **primaQUANT 1STEP PROBE 2x qPCR Master Mix** is subject to strict quality controls. Each lot is tested in a probe-based qPCR with cDNA and DNA input and must conform to our quality control chart.

Enzyme **purity and homogeneity of > 98 %** is validated by Bioanalyzer SDS protein electrophoresis.

The primaQUANT 1STEP PROBE 2x qPCR Master Mix is free of:

- RNA
- DNA
- RNase
- DNase
- Endo- & Exonuclease activity

FURTHER INFORMATION

For more information, please visit our website: www.steinbrenner.de



Steinbrenner
Laborsysteme GmbH

In der Au 17 | 69257 Wiesenbach

 +49 (0) 6223 / 96 73 00

 mail@steinbrenner.de



Further Products

Products that may also interest you

