

Article	Reactions
SL-9802-5X10ML	8x 1.25 mL, 2000 rxn á 25µL

Storage Conditions
Long-Term Storage at -20 °C in the dark
Short-Term Storage at 4 °C in the dark



5x custom qPCR Probe Master Mix - no ROX

DESCRIPTION

Our **primaQUANT PROBE 5x custom qPCR Master Mix** is an optimized ready-to-use mixture for probe-based assays such as Taqman®, Beacons and MGBs. It contains a modified HotStart DNA Polymerase, as well as dNTPS and MgCl₂ - combined in an optimized buffer system that provides fast kinetics and target amplification even for difficult templates.

The **primaQUANT 5x Master Mix** contains all components - you just need to add primers and template DNA/cDNA.

The Master Mix can be used not only for expression analysis but also for genotyping, copy-number analysis and all sorts of probe-based quantitative PCR.



DID YOU KNOW?

- > **primaQUANT PROBE 5x** is compatible with BD MAX™ systems.
- > **primaQUANT PROBE 5x** is especially designed for multiplexing applications.

STANDARD PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Cycling conditions highly depend on the primer, probe, amplicon and input material and thus might need adjustments.
- > However, standard cycling conditions can be applied for the majority of qPCR assays out-of-the box.

3-Step Protocol

Step	Time	Temperature
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 seconds	92 °C - 95 °C
Annealing	5 seconds	60 °C - depending on primer
Extension	10 seconds	72 °C

2-Step Protocol

Step	Time	Temperature
Initial Denaturation	1-3 minutes	92 °C - 95 °C
25- 40 cycles		
Denaturation	5 seconds	92 °C - 95 °C
Annealing/Extension Combined	10 seconds	60 °C - depending on primer

ULTRA-FAST PROTOCOL

BEFORE YOU START

- > After thawing, please **invert the Master Mix tube 6-8 times**.
- > **Do not vortex** the Master Mix to prevent damage of the enzyme.

NOTE

- > Ultra-fast Cycling conditions **highly depend on the ramping rate of your qPCR cycler, primer, probe, amplicon and input material** and thus might need adjustments.
- > Ultra-fast cycling conditions can be applied for the majority of qPCR assays out-of-the box, provided that your primer/probe sets do not show unspecific binding.

3-Step Protocol

Step	Time	Temperature
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1 seconds	92 °C - 95 °C
Annealing	1 seconds	60 °C - depending on primer
Extension	1 second	72 °C

2-Step Protocol

Step	Time	Temperature
Initial Denaturation	60 seconds	92 °C - 95 °C
25- 40 cycles		
Denaturation	1 second	92 °C - 95 °C
Annealing/Extension Combined	1 second	60 °C - depending on primer

RECOMMENDED REACTION MIXTURE PER WELL

Components	25 µL Reaction	10 µL Reaction	Final Concentration
5x primaQUANT Master Mix	5 µL	2 µL	1x
Forward Primer	variable (e.g. 2 µL)	variable (e.g. 1 µL)	100 - 400 nM
Reverse Primer	variable (e.g. 2 µL)	variable (e.g. 1 µL)	100 - 400 nM
Probe	variable (e.g. 2 µL)	variable (e.g. 1 µL)	200 - 600 nM
Template DNA	variable	variable	0.1 - 10 ng/reaction
Sterile Water	adjust to 25 µL	adjust to 10 µL	

NOTE

> For maximum efficiency and specificity annealing temperatures as well as extension time, primer/probe concentration and template concentration need to be optimized.

CALCULATOR TOOL



Please feel free to download our Excel sheet calculator to calculate the necessary volumes:

<http://calculator.steinbrenner-laborsysteme.de>.



<https://www.steinbrenner-laborsysteme.de>

HOW TO VALIDATE A QPCR SETUP


i BEFORE YOU START

- > You can find additional information on how to validate and set up a qPCR in our qPCR Knowledge Center.

Required Controls

- > **DNA Dilution Series**
A DNA dilution series is used to validate the dynamic range, find the optimal DNA input amount and estimate the overall PCR efficiency.
- > **Primer Dilution Series**
High primer amounts can result in unspecific primer binding that limit the fidelity of your qPCR.
- > **No Template Control (NTC)**
A control in which all components except the template are added - this control is used as a negative control and should not show amplification.
- > **No Reverse Transcription (NoRT)**
For this control, reverse transcription is performed without the reverse transcriptase. It helps to identify unwanted amplification from left-over gDNA or unwanted templates.

Recommended Validation Layout



APPLICATIONS

- > Probe-based quantitative PCR
 - TaqMan® Probes
 - Any Dual-Labeled Hydrolysis Probe
 - Molecular Beacons
 - Scorpion Probes
- > DNA Genotyping
- > DNA SNP Analysis
- > RNA and miRNA Expression
- > Multiplexing (up to 4 colors)
- > Transcript Variant Analysis

QUALITY CONTROL

Our **primaQUANT PROBE 5x qPCR Master Mix** undergoes stringent quality controls. Each lot is tested in a probe-based qPCR with cDNA and lambda DNA input.

Enzyme purity and homogeneity of > 98 % is validated using a Bioanalyzer SDS protein electrophoresis.

All Master Mixes are free of detectable endonuclease- & exonuclease activity:

- Incubation of 1 µg of plasmid DNA with 5 U for 4h at 37 °C and 72 °C
- Incubation of 1 µg of a DNA size standard with 5 U for 4h at 37 °C and 72 °C

FURTHER INFORMATION

For more information, please visit our website

<https://www.steinbrenner-laborsysteme.de>



qPCR
KnowledgeCenter



Steinbrenner
Laborsysteme GmbH