

magtivio



MagSi-DNA Animal

*Fast and cost-effective extraction of
genomic DNA from aquaculture samples*



DNA Extraction from Aquaculture Samples

MagSi-DNA Animal allows for fast and cost-effective extraction of DNA from aquaculture samples. We demonstrated the extraction of DNA from shrimp samples for genomics-assisted breeding, from different species of fish fins for genetic testing by aquaculture laboratories, and from Atlantic mackerel for fisheries management. The kit is optimized to extract DNA from different types of aquaculture sample materials with the highest purity and delivering DNA which is suitable for genotypic assays or qPCR based analysis.

Features

- Short (and customizable) protocols, processing at room temperature (after sample lysis)
- Preparation time of 96 samples: 20 min after lysis
- High DNA yield, high molecular weight and excellent purity, suitable for long-time storage
- Suitable for many genomic applications such as SNP genotyping, NGS and qPCR
- Suitable for genomics-assisted breeding, genetic testing and fisheries management
- Can be used on PurePrep, KingFisher™ or other automated DNA purification instruments and compatible with various liquid handling robots (e.g. Hamilton®, TECAN®, oKtopure™)
- Magnetic separators available for manual procedure with microtubes and microplates
- Validated on various species: Atlantic mackerel, rainbow trout, Atlantic cod, Atlantic salmon, Japanese amberjack and black tiger shrimp (check our Application Notes)
- Suitable for processing various sample materials: fish fins, antenna's, pereopods (walking legs), pleopods (swimming legs), and uropods (tails)





Figure 1. DNA concentrations obtained from fins of 4 different fish, measured with the Qubit™ dsDNA BR Assay Kit. The different input amounts (10, 20 and 40 mg) resulted in a DNA yield ranging from 16 to 225 ng/µL. The data presented are mean values (n= 8, ±1 SD). More info in Application Note AN0018.




Figure 2. DNA concentrations obtained from different black tiger shrimp parts, measured with the Qubit™ dsDNA BR Assay Kit. The different shrimp parts resulted in a DNA yield ranging from 21 to 115 ng/µL. The data are presented as mean (n=8, ±1 SD). More info in Application Note AN0019.

Ordering information

Art. No.	Description	Amount
MDKT00150096	MagSi-DNA Animal	96 preps
MDKT00150960	MagSi-DNA Animal	10 x 96 preps
AS00001	PurePrep 96 Nucleic Acid Purification System	1 unit
MDPL00200050	2 mL Deepwell Plate with square wells for KingFisher™/PurePrep 96	50 pcs/pack
MDPL00190060	200 µL square-well Elution Plate for KingFisher™/PurePrep 96	60 pcs/pack
MDPL00210060	96 well Tip-Comb for KingFisher™/PurePrep 96	60 pcs/pack

magtivio B.V.

Daelderweg 9
 6361 HK Nuth | The Netherlands

Tel.: +31 (0)45 208 4810
 Fax: +31 (0)45 208 4817
 www.magtivio.com
 info@magtivio.com

PL0065-001 | © 2022

