



primaQUANT 1-Step RT-qPCR Master Mix with SYBRGreen

Article	Content
SL-9561R_smp	0.5 ml, 50 rxn × 20 µL
SL-9561R-1ML	1 ml, 100 rxn × 20 µL
SL-9561R-5ML	5 ml, 500 rxn × 20 µL
SL-9561R-10ML	10 ml, 1000 rxn × 20 µL
SL-9561R-20ML	20 ml, 2000 rxn × 20 µL



Long-Term Storage
at -20°C in the dark

Short-Term Storage
at 4°C in the dark

DESCRIPTION

Our **primaQUANT 1STEP RT-qPCR Master Mix with SYBRGreen** for probes is an optimized ready-to-use Master Mix for dye-based assays and can be directly used with both DNA and RNA as a starting material.

It contains a modified and proprietary HotStart DNA Polymerase and Reverse Transcriptase, as well as dNTPS, MgCl₂, SYBRGreen and other components in optimized concentrations. It provides fast kinetics, a RNA sensitivity of < 10 fg, target amplification even for difficult templates and multiplexing of more than 6 targets.

The **primaQUANT 1STEP RT-qPCR 2x qPCR Master Mix with SYBRGreen** contains all components - you just need to add primers and template DNA/cDNA or RNA.



DID YOU KNOW?

For qPCR cyclers that do not require ROX as a reference dye, **primaQUANT 1STEP PROBE Master Mix with SYBRGreen** is also available without ROX (SL-9561).



Recommended Reaction Mixture per Well



BEFORE YOU START

- After thawing, please **invert the Master Mix tube 6-8 times**.
- **Do not vortex** the Master Mix to prevent damage of the enzyme.

Component	Stock Concentration	20 μ l Reaction	10 μ l Reaction	Final Concentration
2x primaQUANT 1STEP RT-qPCR Master Mix with SYBRGreen	2x	10 μ l	5 μ l	1x
Forward Primer	4 μ M	1 μ l	0.5 μ l	200 nM (100 - 400 nM recommended)
Reverse Primer	4 μ M	1 μ l	0.5 μ l	200 nM (100 - 400 nM recommended)
Template RNA or DNA/cDNA	-	variable	variable	0.1 - 100 ng / Reaction
Sterile Water	-	adjust to 20 μ l	adjust to 10 μ l	-



NOTE

For maximum efficiency and specificity, adjustments of annealing temperature as well as extension time, primer concentration and template concentration may be needed.

CALCULATOR TOOL

Please feel free to download our Excel sheet calculator to calculate the necessary volumes:
calculator.steinbrenner-laborsysteme.de



qPCR
KnowledgeCenter

Standard Protocol



NOTE

- For the majority of RT-qPCR assays, standard cycling conditions can be applied.
- However, cycling conditions strongly depend on the primer, amplicon and input material and thus some of these factors might need adjustments.

ONE-STEP RT-QPCR WITH ADDITIONAL ANNEALING

	Time	Temperature
Reverse Transcription	10 minutes	50°C
Initial Denaturation	1 - 3 minutes	92°C - 95°C
25 - 40 cycles	Denaturation	5 seconds 92°C - 95°C
	Annealing	5 seconds 60°C depending on primer
	Extension	5 - 10 seconds 72°C
Melting Curve (optional)	-	60°C ramping up to 95°C

ONE-STEP RT-QPCR PROTOCOL WITHOUT ADDITIONAL ANNEALING

	Time	Temperature
Reverse Transcription	10 minutes	50°C
Initial Denaturation	1 - 3 minutes	92°C - 95°C
25 - 40 cycles	Denaturation	5 seconds 92°C - 95°C
	Annealing / Extension combined	5 - 20 seconds 60°C depending on primer
Melting Curve (optional)	-	60°C ramping up to 95°C



Ultra-fast Protocol



NOTE

- Ultra-fast cycling conditions highly depend on the ramping rate of your qPCR cycler, primer, amplicon and input material and thus might need adjustments.
- Ultra-fast cycling conditions can be applied for the majority of qPCR assays out-of-the box, provided that your primer sets do not show unspecific binding.

ONE-STEP RT-QPCR PROTOCOL WITH ADDITIONAL ANNEALING

	Time	Temperature
Reverse Transcription	5 - 10 minutes	50°C
Initial Denaturation	1 minute	92°C - 95°C
25 - 40 cycles	Denaturation	1 second
	Annealing	1 - 5 seconds
	Extension	1 second
Melting Curve (optional)	-	60°C ramping up to 95°C

ONE-STEP RT-QPCR PROTOCOL WITHOUT ADDITIONAL ANNEALING

	Time	Temperature
Reverse Transcription	5 - 10 minutes	50°C
Initial Denaturation	1 minute	92°C - 95°C
25 - 40 cycles	Denaturation	1 second
	Annealing / Extension combined	1 - 5 seconds
Melting Curve (optional)	-	60°C ramping up to 95°C



Applications

Dye-based quantitative PCR with RNA Input

Dye-based quantitative PCR with cDNA/DNA Input

High Resolution Melting Curve

QUALITY CONTROL PROCEDURE

Our **primaQUANT 1STEP PROBE 2x qPCR Master Mix** is subject to strict quality controls. Each lot is tested in a dye-based qPCR with cDNA and DNA input and must conform to our quality control chart.

Enzyme **purity and homogeneity of > 98 %** is validated by Bioanalyzer SDS protein electrophoresis.

The primaQUANT 1STEP PROBE 2x qPCR Master Mix is free of:

- RNA
- DNA
- RNase
- DNase
- Endo- & Exonuclease activity

FURTHER INFORMATION

For more information, please visit our website: www.steinbrenner.de



Steinbrenner
Laborsysteme GmbH

In der Au 17 | 69257 Wiesenbach



+49 (0) 6223 / 96 73 00



mail@steinbrenner.de



Further Products

Products that may also interest you

