


Deep-Well Plates.

primaPLATE Consumables




Steinbrenner
Laborsysteme GmbH

For KingFisher™

High quality tip combs, deep-well plates and elution plates.

- Available as **set** or **separate items**
- RNA / DNA and RNase / DNase-free
- **Individual package sizes** available
- Compatible with:
 - KingFisher™ Flex
 - KingFisher™ 96
 - KingFisher™ Presto
 - Qiagen Biosprint 96
 - Pureprep 96
 - IDEAL 96



Article	Package Size
978-KF-COMB	60 Tip-Combs
978-KF-DEEPWELL	60 Deep-Well Plates
978-KF-ELUTION	60 Elution Plates
978-KF-KIT	individual compositions available



Consumables

Consumables for Liquid Handling.



We also supply **96/384 well plates** for PCR and qPCR, liquid handling reservoirs and filter tips.

All our plates are made from **medical grade ultrapure polypropylene (PP)** and are free from RNA/DNA, bacterial DNA and RNase/DNase.

When it comes to **fast and accurate pipetting** of multiwell plates, labs all over the world trust our manual 96 channel pipettor **Liquidator96®**.



96-Well > 1.1 ml.




This deep-well plate is especially designed for small volumes. It provides round wells with **1.1 ml volume** and a **V-bottom**.

The plate is ideally suited for sample storage, dilution, and mixing. It also perfectly serves as elution plate for magnetic bead applications, e.g. with the TANbead Maelstrom system™.


Article	Max Vol.	Well Type	Bottom	Package Size
978-1100RV	1.1 ml	Round	V-Bottom	60 Plates

96-Well > 2.2 ml.



Our plate with **2.2ml volume** combines square wells with a **U-bottom**.

This plate is ideal for **sample storage** and for liquid handling steps requiring larger volumes such as washing, dilution, or sample preparation. The modified U-bottom is perfect for mixing. The plate is also suited for **magnetic bead** sample processing, e.g. with KingFisher™



Article	Max Vol.	Well Type	Bottom	Package Size
978-2200SU	2.2 ml	Square	U-Bottom	60 Plates



Steinbrenner

Laborsysteme GmbH

In der Au 17
69257 Wiesenbach

+49 (0) 6223 / 96 73 00

✉ mail@steinbrenner.de

🌐 www.steinbrenner.de

Get your
Free Sample & Support.

🌐 » www.steinbrenner.de



Manfred Lux



Uli Christ