

Data sheet

Colli-Pee® UAS™ FV-5010-RUO



Product description

Product Name:	Colli-Pee® UAS™ FV-5010-RUO
Reference Number:	N00345_RUO
Components:	Housing, floater, IFU, tube cap, collector tube, and UAS preservative
Order Quantity:	25 or 100

Product Characteristics

Material:	Polypropylene (PP) and High Density Polyethylene (HDPE)
Packaging:	See-through pouch constructed of medical grade paper/multiply Polyethylene Teraphtalate / Polypropylene (PET/PP) plastic laminate
Product Dimensions:	240 x 135 x 40* mm
Product Weight	± 38 g
Tube Volume:	13 mL
Collected Volume:	9-12 mL
Preservative Volume:	3 mL
Collector Tube - Tube Cap Interface:	Screw thread
Housing - Collector Tube Interface:	Push-fit
Shelf-life:	1 year
Storage Temperature:	15-30°C
Sterility:	Non-sterile



*The height of the pouched device can be squeezed to 30 mm when being included in a postal kit.