

# Neptune<sup>®</sup> Filter Tips.

*Universal Filtered Pipette Tips*



**Steinbrenner**  
Laborsysteme GmbH

# Neptune<sup>®</sup> Universal Filter Tips.

Our Neptune<sup>®</sup> filter tips meet highest quality expectations and guarantee best pipetting results. They come in clear and sturdy racks with hinged lids and are certified free of RNA/DNA, RNase/DNase and endotoxins.

- **FlexFit technology** – fits most common pipettes
- **S3 low retention** polymer – no coating
- **Extra-long** 10 µl tip
- **Extra-long** 1250 µl tip

Article	Max. Volume	Package Size
976-010	10 µl	5 x 960 tips (5 x 10 racks)
976-010XL	10 µl – extra long	
976-020	20 µl	
976-100	100 µl	
976-200	200 µl	
976-300	300 µl	4 x 768 tips (4 x 8 racks)
976-1000	1000 µl	
976-1250	1250 µl – extra long	



Consumables

## S3 Polymer.

The excellent **sample-saving surface (S3)** property of our Neptune® filter tips guarantees most accurate pipetting results. The superior low-retention characteristics are achieved without additional treatment or surfactant agents.

Neptune® filter tips are manufactured in state-of-the-art cleanrooms to meet **highest quality standards.**

## Optimized Tip Design.

The **extra-long versions** of our 10 µl and 1250 µl Neptune® filter tips are perfect for the handling of samples stored in tubes and plates. With its outstanding length of 9.75 cm, the 1250 µl tip allows to reach the bottom of almost every tube.

Their **slim shape** minimizes the danger of liquid overflow which makes these tips especially suited for narrow tubes and plates.

**NEPTUNE**  
Tools for Life Sciences

# Universal Fit.

Our Neptune® filter tips with FlexFit sealing technology are designed to **fit onto most pipette types** on the market.

For evaluation with your pipettes, please contact us for a free sample.

976-010



976-010XL



976-020



976-100



976-200



976-300



976-1000



976-1250



# Increased Sustainability.

Our Neptune® product line has been completely redesigned to address **environmental aspects**.

The improved rack **saves 30% plastic** and is made from clear material to facilitate recycling.





# Steinbrenner

Laborsysteme GmbH

In der Au 17  
69257 Wiesenbach

+49 (0) 6223 / 96 73 00

mail@steinbrenner.de

www.steinbrenner.de



Get your  
**Free Sample & Support.**

 » [www.steinbrenner.de](http://www.steinbrenner.de)



Dr. Jan Winter



Manfred Lux