


# Plates for PCR & qPCR.

*primaPLATE Consumables*



**Steinbrenner**  
Laborsysteme GmbH

# Plate Seals.

High-quality adhesive seals and heat seals for your PCR and qPCR applications.



Article	Seal Characteristics
SL-PS-PCR	Adhesive seal, clear, DMSO-resistant
SL-PS-PCR-F	Adhesive seal, aluminium foil, strong
SL-PS-PCR-C	Adhesive seal, clear, strong
SL-PS-qPCR	Adhesive seal, optically clear, self-adhesive
SL-PS-qPCR2	Adhesive seal, optically clear, pressure-sensitive
SL-PS-qPCR3	Adhesive seal, optically clear, pressure-sensitive, wide temperature range
SL-PS-HS1	Heat seal, piercable, resealable
SL-PS-HS6	Heat seal, heavy-duty

More seals can be found on our website [www.steinbrenner.de](http://www.steinbrenner.de).



Consumables

 **prima**PLATE

# High-Quality PCR & qPCR plates.

- Made of medical polypropylene
- Certified free of RNA/DNA, RNases/DNases and endotoxins
- Versions for most PCR/qPCR cyclers available

 **prima**PLATE



# Premium-Quality PCR & qPCR plates.

- Highest precision for demanding tasks
- Perfect for automation due to reinforced frame
- Two-component design guarantees high precision and rigidity and minimal evaporation losses
- Wells made of medical polypropylene
- Fiberglass reinforced frame
- Certified free of RNA/DNA, RNase/DNase
- Certified free of heavy metals and PCR inhibitors
- Certified animal-free
- Version for most PCR/qPCR cyclers available

# Compatible with various cyclers.

Manufacturer	Type	Plate
Applied Biosystems / Thermo Fischer	Quantstudio 3/5/7/12	SL-PP96-1, SL-PP96-2, SL-PP96R-2CC*
	Quantstudio 3/5/7/12 FAST Block	SL-PP96-1L, SL-PP96-2L, SL-PP96-3L
	Quantstudio 3/5/7/12 384-Well	SL-PP384-SC, SL-PP384R-BC*, SL-PP384R-WC*
Analytik Jena	qTower3 96-Well	SL-PP96-1*, SL-PP96-1L, SL-PP96-2, SL-PP96R-2-WW*, SL-PP96-2L, SL-PP96-LC, SL-PP96R-4-CC*, SL-PP96-4L
	qTower3 384-Well	SL-PP384-SC, SL-PP384R-BW*, SL-PP384R-WW*
BioRad	CFX96 Touch	SL-PP96-1L, SL-PP96-2L, SL-PP96-4L, SL-PP96R-4L-WW*
	CFX384 Touch	SL-PP96-SC, SL-PP384-DC, SL-PP384R-BW*
	CFX Opus 96	SL-PP96-2L, SL-PP96-4L, SL-PP96R-4L-WW*
	CFX Opus 384	SL-PP96-SC, SL-PP384-DC, SL-PP384R-BW*
Agilent Technologies	AriaMx	SL-PP96-1L, SL-PP96-4L, SL-PP96R-4L-WW*
Eppendorf	Mastercycler ep	SL-PP96-1*, SL-PP96-1L, SL-PP96-2, SL-PP96R-2-CC*, SL-PP96R-4-CW*, SL-PP96-4L, SL-PP96R-4L-WW*
Roche	Lightcycler 480	SL-PP96-LC, SL-PP96R-LC-WW*

\* available in various frame and well colours

All variants and details can be found on our website [www.steinbrenner.de/plates](http://www.steinbrenner.de/plates)



# Steinbrenner

Laborsysteme GmbH

In der Au 17  
69257 Wiesenbach

+49 (0) 6223 / 96 73 00

✉ [mail@steinbrenner.de](mailto:mail@steinbrenner.de)

🌐 [www.steinbrenner.de](http://www.steinbrenner.de)

Get your  
**Free Sample & Support.**

🌐 » [www.steinbrenner.de](http://www.steinbrenner.de)



Manfred Lux



Uli Christ